





# QUICK, EFFICIENT SOLUTIONS FOR

Configure it to alert you when a patient/resident:

- Sits up in their bed
- Leaves their bed Leaves their room
- Leaves their chair (chair sensor plate/pad can be purchased separately)
- Set a time delay for chair/bed plate (Time delay 0, 1 sec, 2 secs, 3 secs, 10 secs, 10 mins, 20 mins) for more independent residents that can perform certain tasks unaided (like visiting their en-suite)

A unique, easy to use unit that's reliable, portable and cost-effective.






PLATE OR PAD HERE

**PLUG AND PLAY** Will work with your Nurse Call system or your money back

### UNIQUE, FLEXIBLE, PORTABLE

## ONE OF A KIND

Dementia, frailty and confusion are amongst the most common reasons for resident injury, to name but a few. Different residents have specific needs which currently necessitates the purchase of a wide range of peripherals. Live-Link's Movement Sensor makes carers lives easier by allowing you to use one device for a multitude of purposes in one, robust, easy to use, reliable unit. It's one of a kind. The Live-Link Movement Sensor can be set to register a call system alert under a range of different needs.

be to monitor residents at risk of falling. Unlike a pressure mat, our sensor can't be stepped over or avoided. We feel that the Live-Link Movement Sensor could put an end to these risks.

Over the years we have listened to feedback from customers about how hard it seems to

**FREE TRIAL** 



# SENSOR

WHEN A RESIDENT ATTEMPTS TO LEAVE THEIR BED WHEN A RESIDENT HAS LEFT THEIR BED WHEN A RESIDENT EXITS THEIR ROOM

\* Based on 4 Duracell Industrial AA batteries and up to 10 activations per day

### **PARTNER PRODUCT**

**LIVE-LINK** 

steel chair/bed plate can be connected to your Movement Sensor. The Live-Link Sensor Plate is ideal for monitoring somebody either sitting in their chair or as a bed sensor. Three year Warranty

The perfect accessory to aid fall prevention, the Live-Link

- Replaceable Hygienic Plastic Covers Stainless Steel

• Links into Movement Sensor

## SENSOR PLATE

LIVE-LINK

**EFFORTLESS** OPERATION Live-Link's Movement Sensor is fully portable and can be moved very easily from one location to another and can be linked into the majority of nurse call systems. The unit

comes pre-programmed with

4 x AA Batteries fitted and

can be installed and

operational within a

few minutes.

Low-battery indicator: The red LED on the Movement Sensor will be permanently lit. Repeated calls are sent to the nurse call unit, which it is linked into, to alert staff that new batteries are required. All four batteries need to be replaced at the same time. Please note this piece of equipment is designed to compliment the care of the patient and not replace good nursing care. Live-Link has a policy of continuous product development and reserves the right to change any prices and specifications given in this brochure.

Movement Sensor is powered by 4 x AA Batteries. Battery Life: 12 months to 3

years depending on usage and battery type

Other products from Live-link

\*\* Movement Sensor Only (does not cover faults arising out of neglect or careless handling)



**FREEPHONE:** 

08080 687732

E: livelink@nursecallalarms.co.uk

CHAIR PADS

WARRANTY ON

**PARTS AND REPAIR**\*\* www.nursecallalarms.co.uk

FLOOR MATS